




## iGEM Pasteur Paris 2016 Protocols



A CUSTOMIZABLE  
BIODETECTION SYSTEM

# DEPHOSPHORYLATION

**Aim:** To reduce vector self-ligation and favor that of vector to insert instead.

The presence of a free phosphate group at one end of the linear vector DNA after a restriction enzyme digestion will lead to self ligation, because it is intra-molecular.

### Materials:

- Molecular biology equipment: 37 °C and 65 °C water baths, pipette set, ice bucket
- (Recombinant Shrimp alkaline phosphatase) rSAP Reaction Buffer (10X)
- rSAP
- DNA
- Nuclease-free water

### Protocol:

1. Digest vector with the appropriate [restriction enzymes](#). Enzymes that leave non-compatible ends are ideal as they prevent vector self-ligation as well.
2. Dephosphorylation is necessary to prevent self-ligation. NEB offers three products for dephosphorylation of DNA:
  - rSAP ([NEB #M0371](#)) is heat-inactivable phosphatases -> It works in all NEBuffers, but rSAP requires that the reaction be supplemented with  $Zn^{2+}$

### **Dephosphorylation of 5' ends of DNA using rSAP:**

- o rSAP Reaction Buffer (10X) : 2  $\mu$ l
- o DNA  $\geq$  1 pmol of DNA ends (about 1  $\mu$ l of 3 kb plasmid at 100 ng/  $\mu$ l)
- o 1  $\mu$ l rSAP (1 unit/  $\mu$ l)
- o Nuclease-free Water : to 20  $\mu$ l
- o Incubation : 37°C for 30 minutes
- o Heat Inactivation : 65°C for 5 minutes

Note: Scale larger reaction volumes proportionally

---